

BioSentry® Chlorhexidine Solution

CONTAINS 2 % CHLORHEXIDINE GLUCONATE

ACTIVE INGREDIENT:

Chlorhexidine Gluconate.....2%

EPA Est. No. 61282-WI-01

**FOR EXTERNAL USE ONLY
KEEP OUT OF REACH OF CHILDREN**

INDICATIONS: A topical aqueous cleaning solution.

CAUTION: Avoid contact with eyes and the mucous membranes. This product is not to be used in ears. If contact is made, flush immediately and thoroughly with clean water.

WARNINGS: For external use only.

DOSAGE AND ADMINISTRATION:

Dilute one (1) once (2 tablespoons) of Chlorhexidine solution per gallon of clean water. Rinse area to be disinfected with ample amount of Chlorhexidine solution. Wipe away the excess and pat dry with a sterile gauze or sponge.

Store at controlled room temperature between 15°-30°C (50°-86°F). Keep container tightly closed when not in use.

**Manufactured By:
HACCO, Inc.
110 Hopkins Drive
Randolph, WI 53956 U.S.A.**

BioSentry® is a Registered Trademark of BioSentry, Inc.; used under Exclusive License by Neogen Corporation.

Code: NA.RR.

Item No.: D190450
Net Contents: 1 Gallon (3.785 Liters)
UPC: 095242 02400(6)
Format: UPC-A
Label: LA 220010

Item No.: D190450
Net Contents: 4 x 1 Gallon
UPC:
Format: Interleaved 2-of-5
Label: Case Label