

BioSentry® LongLife 250S™

SUPREME SPECTRUM VIRUCIDAL DISINFECTANT

For Extended Protection against Avian Diseases

COMPOSITION:

An active synergistic blend of organic acids, organic biocides and surfactants.

Reg. No.: To Be Assigned
EPA Est. No. 61282-WI-1

**KEEP OUT OF REACH OF CHILDREN
DANGER**

APPLICATION

Heavy duty disinfectant, ideal for high challenge environments. For farm disinfection of buildings, equipment, pathways and aprons, including foot-dips and wheel-dips, as part of a continuous biosecurity program. Effective against all families of virus as well as mycoplasma, bacteria, fungi (including *Aspergillus fumigatus*), yeasts and moulds.

DIRECTIONS FOR USE

Use at a dilution of 1:90 – 1:400 (110 ml – 25 ml in 10 litres of water) depending on the level of organic challenge.

Application Rate: 300 ml disinfectant solution per 1 m² of surface area.

Precautions: Ensure adequate ventilation to prevent vapour build-up. Store tightly capped in a cool place. Keep in original container. Do not store or produce near dairy products. Avoid contact with strong oxidizing agents and bases.

Disposal of Waste: Do not allow product to enter surface waters and ground waters. Dispose of any waste product in accordance with the safety data sheet. This material must be disposed of as special waste. Always consider local restrictions.

BioSentry® LongLife 250S™ Contains: Tar acids, polyalkylphenol fraction; Benzenesulphonic acid, mono-C₁₀₋₁₄-alkyl derivs.; Chlorinated xylenols fraction.

CORROSIVE

- Harmful in contact with skin and if swallowed
- Causes burns
- May cause sensitization by skin contact
- Harmful to aquatic organisms
- Keep away from food, drink and animal feedingstuffs
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
- Wear suitable protective clothing, gloves and eye/face protection
- In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible)

- Avoid release to the environment. Refer to special instructions/safety data sheets

**Medical 1-651-523-0318 (International)
? 1-800-621-8829**

**Manufactured For:
HACCO, Inc.
110 Hopkins Drive
Randolph, WI 53956 U.S.A.**

**Imported By:
Oscar Francois Limited
Alicia House
135 Duke Street
Port of Spain
Trinidad and Tobago**

BioSentry® is a Registered Trademark of Neogen Corporation.
LongLife 250S™ is a Trademark of Neogen Corporation.

Item No.:
Net Contents: 3.785 Liters (1 Gallon)
UPC:
Format:

Label Code: NA.NO.